

# ARION

## PTA Balloon Catheter

Balloon Catheter for treatment of Peripheral Arterial Stenosis



### EFFORTLESS CROSSING WITH GOOD SUPPORT

- Over the wire (OTW) constructions, with 0.035" guidewire for better support
- Surmodics™ hydrophilic coating offers excellent gliding
- Tip profile of 0.045" with a tapered profile for easy entry

### RANGE OF BALLOON OPTIONS

- Wide range of lengths for treatment of long, diffuse lesions
- Available in 4-10 mm diameter and 20-100 mm in length
- Can be effortlessly used in both above-the-knee and below-the-knee applications

### SUPERIOR BALLOON CONFORMANCE


- Specially selected nylon polymer for superior compliance
- Nominal pressure of 7 atm and Burst pressure of 15 atm gives with a good anatomical compliance

### Arion Balloon Features

Balloon Dia (mm)	Balloon Length (mm)	Catheter Length (cm)	Introducer Sheath (Fr)	Nominal Pressure	Rated Burst Pressure
4.0	20, 40, 60, 80, 100	80 cm	5 Fr	8 atm	14 atm
5.0					
6.0					
7.0		135 cm	6 Fr	6 atm	12 atm
8.0					
9.0					
10.0					

# ARION

## PTA Balloon Catheter



### ORDERING MATRIX

Balloon (mm)	Catheter Length (cm)	BALLOON WORKING LENGTH (MM)				
		20	40	60	80	100
4	80	OTW0200400P14-080	OTW0200800P14-080	OTW0200120P14-080	OTW0200150P14-080	OTW0200150P14-080
5		OTW0500200P35-080	OTW0500400P35-080	OTW0500600P35-080	OTW0500800P35-080	OTW0500100P35-080
6		OTW0600200P35-080	OTW0600400P35-080	OTW0600600P35-080	OTW0600800P35-080	OTW0500100P35-080
7		OTW0700200P35-080	OTW0700400P35-080	OTW0700600P35-080	OTW0700800P35-080	OTW0500100P35-080
8		OTW0800200P35-080	OTW0800400P35-080	OTW0800600P35-080	OTW0800800P35-080	OTW0500100P35-080
9		OTW0900200P35-080	OTW0900400P35-080	OTW0900600P35-080	OTW0900800P35-080	OTW0500100P35-080
10		OTW1000200P35-080	OTW1000400P35-080	OTW1000600P35-080	OTW1000800P35-080	OTW0500100P35-080
4	135	OTW0400200P35-135	OTW0400400P35-135	OTW0400600P35-135	OTW0400800P35-135	OTW0400100P35-135
5		OTW0500200P35-135	OTW0500400P35-135	OTW0500600P35-135	OTW0500800P35-135	OTW0500100P35-135
6		OTW0500200P35-135	OTW0500400P35-135	OTW0500600P35-135	OTW0500800P35-135	OTW0500100P35-135
7		OTW0500200P35-135	OTW0500400P35-135	OTW0500600P35-135	OTW0500800P35-135	OTW0500100P35-135
8		OTW0500200P35-135	OTW0500400P35-135	OTW0500600P35-135	OTW0500800P35-135	OTW0500100P35-135
9		OTW0500200P35-135	OTW0500400P35-135	OTW0500600P35-135	OTW0500800P35-135	OTW0500100P35-135
10		OTW0500200P35-135	OTW0500400P35-135	OTW0500600P35-135	OTW0500800P35-135	OTW0500100P35-135

(\* sizes are available only on demand)

Distributed by

**bioradmedisys**  
science for people

bioradmedisys.com  
Biorad Medisys Pvt. Ltd.

Follow us on    

Corporate Office

Survey No.: 48/3 & 48/7, Pashan Sus Rd, Sus Village, Taluka Mulshi, Pune - 411021. MH, India.

# ARION

---

PTA Balloon Catheter