





VESPREP

PTA Balloon Catheter

Crafted from Physician Insights

VESPREP

PTA Balloon Catheter

Features

RBP upto

20 atm

for optimal vessel prep

Catheter length

60 cm

use with single length wire

Hydrophilic Coating


upto **30** cm

for smooth tracking in the vessel

Balloon length

40 mm

for optimal stenosis dilation



Ordering Information

VESPREP

PTA Balloon Dilation Catheter, GW 0.035"

Product Code	Balloon Diameter (mm)	Balloon Length (mm)	Catheter Length (cm)
OW3505040S	5	40	60
OW3506040S	6	40	60
OW3507040S	7	40	60

Compliance Chart

	2	4	6	8	10	12	14	16	18	20
5	4.2	4.6	5	5.1	5.3	5.5	5.6	5.9	6.1	6.4
6	4.8	5.5	5.8	6.1	6.3	6.5	6.7	6.9	7.2	7.4
7	6.1	6.7	7.1	7.4	7.6	8.0	8.2	8.5	8.9	9.3
	NP*									RBP**

*NP: Nominal Pressure

**RBP: Rated Burst Pressure

bioradmedisys[™]
science for people

www.bioradmedisys.com
www.neurovasive.com

Biorad Medisys Pvt. Ltd.

Survey No.: 48/3 & 48/7, Pashan Sus Rd, Sus Village;
Taluka Mulshi, Pune - 411021, Maharashtra, India.
T: +91-020-3091 2000 F: +91-020-3091 2050

NEUROVASIVE
Extension of life

BRP/CAT/NV/VSP/01/00;10-24